



DEFIBRILLATOR CHECKLIST


Every year there are about 30,000 cardiac arrests in the UK that happen outside of a hospital environment. For every second that defibrillation is delayed, a casualty's chances of survival falls by about 10%. We've put together a few tips to help you ensure your defibrillator is always ready for use.

- 1** Set up a regular inspection procedure for your defibrillator.
- 2** Record your inspections. (We've put together a simple sheet to help you).
- 3** Make a note of the expiration date of pads and batteries.
- 4** Reorder batteries and pads before they expire.
- 5** Follow the guidance in your product manual when cleaning your defibrillator

Defibrillator ready to use? Y/N	Date	Initial

Batteries expiration date

Pads expiration date



† Calls cost 7p per minute plus your company's access charge. Excludes weekends Monday-Friday 8.30am to 5.30pm